



# UnityMol – SweetUnityMol

## User Manual

### Installing the software

The software described below has been developed based on revision 676 of UnityMol and version 4.5.2f1 of Unity3D. All source code was implemented using C# and Cg languages built into Unity 3D and is available along with executables for Mac, Windows and Linux platforms on the sourceforge project website.

<http://unitymol.sourceforge.net/sweet/>.

Documentation, input files, and series of illustrations can be found at <http://glycopedia.eu/SweetUnityMol/>.

### Windows version (date: 2014 10 02)

#### 32 bits

umol-win32-20141002\_Data  
Run\_SweetUmol\_32.bat  
umol-win32-20141002.exe

#### 64 bits

umol-win64-20141002\_Data  
Run\_SweetUmol\_64.bat  
umol-win64-20141002.exe


### Mac version (date: 2014 10 02)

umol-macosx-universal-20141002.app

### Linux version (date: 2014 10 02)

umol-linux-20141002\_Data  
umol-linux-20141002.x86  
umol-linux-20141002.x86\_64

## Running Windows version 32 and 64 bits



Our shader need OpenGL. Use the file **Run\_SweetUmol\_32.bat** or **Run\_SweetUmol\_64.bat** depending your version (the line has to be something like "umol-winXX-20141002.exe – force-opengl").

**DO NOT USE THE FULL SCREEN MODE** (check "**Windowed**" in the "UnityMol Configuration" menu (this menu is the menu where you choose your resolution and click on "Play!").

## Description of the coordinate input file (see Annex)


## Ring Color Coding (see Annex)

## Main Menu



Some of these options offering advanced visualization and manipulations (haptic arm) are being implemented and are not functional at this time: October 2010

### 1 – File (PDB Loading)



1 – Open a file from hard drive (you can't change hard drive yet).

**UnityMol and the pdb file have to be on the same drive but on different directories).**

2 – Fetch a pdb from the pdb database (if "Proxy Server" and "Proxy Port" are not empty, please remove all fields otherwise the program will crash)

By default the visualization of UnityMol is the Particle mode. This mode uses the particle system of Unity3D to display the molecule. This may not be the most appropriate visualization. Changing to SmoothHyperballs mode offers a better visualization.


### 2 – Atom



Different types of representations are available by **Atom** or by **Bond**.

The preselected visualization called **Hyperball** (the most interesting one).

Hyperball Visualization



The **Renderer** and **Panels** menus are made to change the **color** and the **texture** of either the entire molecule, or a selection.

**Hide:** Hide the molecule (This option is activated when some visualization

option (like secondary structures) is activated. Uncheck this box to display both the secondary structure and the molecule.


**Radius:** Change the radius of spheres.

**LOD mode :** when you are in Hyperballs mode and this option is checked, when you will change the orientation of the camera the molecule will go back to the *particle* mode and then recreate the Hyperballs. This is to be used for large macromolecular systems for which the Hyperballs representation reduces frame rate too much. FPS: Frame Per Second.

**Automove** (or SPACE key): the camera will turn around the molecule.

**Lock Camera:** The camera will not be able to move when this option is checked.

**Brightness:** Adjust the brightness of Hyperballs.



### [A – HyperBall Style menu \(3\).](#)

**Shrink:** To adjust the shape of the bond.

**Scale:** To adjust the scale of the bond.

**DFactor:** only for Biological Networks.

**Metaphor:** pre-configured style (CPK, Licorice, VdW, Smooth Hyperballs)

**Interactive mode:** use spring forces to make an interactive molecule (An atom can be dragged and he will go back to his original position).


**Drag/Spring:** spring parameters for interactive mode.

**Velocity Colors:** colors each atoms in interactive mode with the velocity of each one.


**Measure dist:** activate and click on 2 atoms to display their distance.

**Dist. Cueing:** adjust the light according to the distance of the camera.

### [Exemple Ball & Sticks, Licorice, VdW, Smooth Hyperballs](#)



### B - Rendered Menu (4)



The whole molecule is selected by default.  
The selection can be changed by:

- Atom Type (A)
- Atom name (B)
- Residue name (C)
- By chain (D)

You can combine your selection.

You have pre-configured coloration on the **Panels** menu (E), and you can change the color of your selection by clicking on the Color Button (F).

A texture can be applied to the selection (G) on your selection. Several textures are available on other pages (you can switch pages by clicking "<<" and ">>").

### C – Panels Menu (3)



This menu contains pre-configured colors. You can switch by clicking on the buttons.

#### Color Panels:


- **All white:** make the molecule all white.
- **Goodsell:** make the color softer.
- **WaterColor:** color carbons in blue.
- **Pastel:** default colors.
- **CPK:** color with the CPK color code.
- **Basic:** color carbon in green.
- **IUPAC:** color with the IUPAC color code.
- **Chains:** color by chain
- **Sugar:** color by sugar Type (the molecule must contain sugar).
- **ADN/ARN:** coloration for DNA and RNA.

#### Texture Panels:

- **Default:** default texture

- **CPK:** CPK coloration with adapted texture.
- **Acid-Basic:** coloration and texture according to the acid-basic status of amino acids.

### 3 – Secondary Structures



Activate the Secondary Structure visualization

Adjust parameters for the secondary structure visualization.

Color by secondary structure type.

Color by chain type

Don't forget to **click on this button to apply changes**.

Trace the shape of the proteins by a spline passing by all carbon alpha. Smoothness is used to adjust the smoothness of the trace.

Like the C-alpha trace, but adjust the radius and the color of each sphere with the B factor value.

Adjust the scale of the B-factor used.

### 4 – Surface



**T:** Density threshold used to generate the isosurface

**Generate:** generate the surface.


**BFactor:** adjust the shape of the surface according to the bfactor value.

**Volumetric:** show the volumetric space of the protein with a particle system.

**Toggle surface:** hide surface representation (when you change a parameter, you have to hide the surface and press the *Generate* button again).

**HetAtoms / Sugars :** show surface for HetAtoms or/and Sugars molecules.

Once you press the *Generate* button, this window will appear.



**Color:** change the color of the surface (external face).

**Inside color:** change the inside color.

**Use atom color:** Color the surface according to the type of atom under each portion of it.

**Use chain color:** use the color of the chains.

**Hydrophobic scale:** Color the surface according to different hydrophobic scale of residues.

**Use properties color:** Color the surface according to the properties of amino acids (basic, acids etc...)

**Use BFactor color:** Color the surface according to the Bfactor value of the surface atoms.

**Texture:** apply a texture on the surface.

**Static cut:** cut the surface (you can control the cutting by a pad that will appear in the upper right corner).

**Mobile cut:** cut the surface and the cut will be the same when you change the camera's orientation.

**Brightness:** adjust the brightness of the surface

**Color weight:** adjust the color concentration on the surface.

## 5 – Electrostat.

To use this menu, you have to have in the same folder as the PDB file (and the same name) a file “.dx” and “.apf” generate with softwares like APBS (Software for evaluating the electrostatic properties of biomolecular systems).



**Transparency:** Add transparency to the mesh generated when you load the negative and positive electrostatic visualization.

**T:** Threshold used to generate the electrostatic isosurface.

**Load Neg/Pos:** load electrostatic visualization.

**Toggle Neg/Pos:** show or hide the electrostatic visualization.

**Volumetric Fields:** volume rendering of the electrostatic fields.

**Field Lines:** show the Field Lines visualization (examine the local intensity of electric fields).

## 6 – Display



This menu is made to **take screenshot**, change the background and add special effects.

**Screenshot** (or P key): take a screenshot (prefer the P key or you will get all menus with the capture)

**ScreenShot Sequence**: doesn't work yet.

**Background**: add a picture in the background (switch between yes/no)

**White/Gray/Black**: quick selection of background color

**Background Color**: choose your background color with a color panel.

**Effects**: add special effect.

**Infos**: show/hide FPS info and Atom/Bond number.

List of all visual effects:

- **SSAO**: Screen Space Ambient Occlusion. Works on spheres, cubes, surfaces and secondary structures
- **BLUR**: add a blur effect when you move the camera.
- **NOISE**: add noise on the screen
- **BLUR2**: constant blur.
- **DOF**: focus on a point and make fuzzy all atoms far from this point. To use you have to be in "Sphere" mode (Atoms -> Atom Style -> Sphere), and click on an atom. Works better for bigger molecules.
- **CREASE**: doesn't work yet.
- **EDGE**: the edge of every atom is black.
- **VORTEX**: twist the molecule around the center.
- **GRAYS**: make the molecule black and white.
- **TWIRL**: similar to VORTEX
- **SEPIA**: use a SEPIA coloration.
- **GLOW**: make all atoms shiny.

## 7 – Advanced

This menu offers advanced options (some are being implemented)



**GUI Scale**: change the size of the GUI.

**Ortho/Persp**: doesn't work yet.

**Best Textures**: don't use anymore.

**Depth Cueing**: doesn't work yet.

**Volumetric Depth Cueing**: Doesn't work yet.


**Ambient Occlusion**: darkens the densest parts of the molecule



to improve depth perception.

## 8 – Guided Navigation

This menu is used to make a “guided navigation”, but for now it works only with the pdb 3EI0 (GLIC) as part of work in progress.



**Symmetry origin:** X,Y,Z coordinates of the center of symmetry.

**Symmetry direction:** vector of the symmetry direction.

**Send:** activate the guided navigation mode.

## 9 – Sugar



This menu is made to deal with the specific features of sugar visualization (principally made with the POLYS 2.0 software. But the visualization mode called *Ring Blending* works with all molecules.

**Enable Ring Blending:** enable the filling of all rings (aromatic, sugar, or other cycle) with a semi transparent color.

**Enable SugarRibbons:** enable the SugarRibbons visualization (to transform a sugar into a “schematic” representation like secondary structures, but adapted for sugars).

**Hide Hydrogens:** hide hydrogens atoms

**Sugar:** hide sugar atoms

**Non Sugar:** hide non sugar atoms.

### Tune Menu

**Show Oxygens:** to display the intracyclic oxygen atom as a sphere.

**Sugar Only:** apply the sugar Ribbons visualization only on “sugar residues” (will detect 3letters sugar code in the PDB file. List of recognize 3 letters word bellow).

**Change Coloration:** open the Color Tune Menu to change color of : sphere, rings and outer ring.

**Oxygen Sphere size :** To adjust the size of the radius of the sphere.

**Ribbons Thickness:** to adjust the thickness of the ribbon.


**Inner Ring Thickness:** to adjust the thickness of the ring.

**Outer Ring Thickness:** to adjust the thickness of the ring.

**Pyranose:** to adjust the thickness of the C1,C4 bond (for pyranose).


**Pyranose:** to adjust the thickness of other bonds (for pyranose).

**Furanose:** to adjust the thickness of bonds for Furanose.



**Do not forget to click on Apply changes to apply the selected new parameters.**

## Color Tune Menu



### For Oxygen Sphere, Outer Rings and Bond:

**Sugar:** use the sugar color to color Oxygen/Outer Rings or Bonds.

**Chain:** use the chain color to color Oxygen/Outer Rings or Bonds

**Custom** color: open a color picker windows to choose a custom color for Oxygen/Outer Rings or Bonds.

You can make the color **darker** or **lighter** with the **Light Color Factor** slider.

*Apply Changes button not needed for this feature.*

## 10 – VRPN

This menu is used to connect a **haptic arms**, but you need a VRPN server installed on your computer. Work in progress.

## 11 – MDDriver

This menu is used to connect a molecular simulation in real-time. Work in progress.

## 12 – Reset

**Reset the view.**

## Binding Keys

LEFT Arrow – A	Move LEFT
RIGHT Arrow – D	Move RIGHT
DOWN arrow - S	Move DOWN
UP arrow – W	Move UP
Q	Rotation DOWN
E	Rotation UP
X	Rotation LEFT
N	Rotation Right
B	ZOOM +
P	ZOOM –
P	Screenshot (only on local pdb)
BACKSPACE	Hide GUI
SPACE	Automove
R	Reset to center
C	Center to one or several atoms (select an atom with the selection mode).

### *Guided Navigation (only for 3EIO)*

U	Constrained Navigation UP (panoramic mode)
J	Constrained Navigation DOWN (panoramic mode)
I	panoramic mode
Numpad + or B	Spreading chains
Numpad – or V	Narrowing chains
Numpad * or T	Reset chains

































## Format of the PDB File

```

ATOM      1  C1  FUC      1      2.387 -19.488 -11.522  0.00  0.00  C1
ATOM      2  C2  FUC      1      3.248 -19.891 -12.734  0.00  0.00  C2
ATOM      3  C3  FUC      1      4.696 -19.417 -12.550  0.00  0.00  C3
ATOM      4  C4  FUC      1      5.229 -19.916 -11.196  0.00  0.00  C4
ATOM      5  C5  FUC      1      4.278 -19.476 -10.072  0.00  0.00  C5
ATOM      6  C6  FUC      1      4.703 -19.967  -8.682  0.00  0.00  C6
ATOM      7  OR  FUC      1      2.950 -19.963 -10.295  0.00  0.00  O7
ATOM      8  OW  FUC      1      2.124 -18.089 -11.437  0.00  0.00  O8
ATOM      9  O2H FUC      1      2.615 -19.391 -13.909  0.00  0.00  OH2
ATOM     10  O3H FUC      1      5.479 -19.850 -13.659  0.00  0.00  OH3
ATOM     11  O4H FUC      1      5.336 -21.339 -11.160  0.00  0.00  OH4
ATOM     12  H5  FUC      1      4.272 -18.370 -10.039  0.00  0.00  H5
ATOM     13  H1  FUC      1      1.402 -19.965 -11.595  0.00  0.00  H1
ATOM     14  H2  FUC      1      3.255 -20.998 -12.784  0.00  0.00  H2
.....
ATOM     23  C1  GAL      2      1.456 -16.523  -9.755  0.00  0.00  C1
ATOM     24  C2  GAL      2      1.018 -17.688 -10.638  0.00  0.00  C2
ATOM     25  C3  GAL      2     -0.154 -17.316 -11.531  0.00  0.00  C3
ATOM     26  C4  GAL      2     -1.262 -16.689 -10.717  0.00  0.00  C4
ATOM     27  C5  GAL      2     -0.709 -15.522  -9.900  0.00  0.00  C5
ATOM     28  C6  GAL      2     -1.722 -14.822  -9.021  0.00  0.00  C6
.....
ATOM     44  C1  XYL      3      2.801 -16.388  -6.544  0.00  0.00  C1
ATOM     45  C2  XYL      3      3.027 -16.035  -8.001  0.00  0.00  C2
ATOM     46  C3  XYL      3      4.524 -16.021  -8.286  0.00  0.00  C3
ATOM     47  C4  XYL      3      5.231 -15.028  -7.349  0.00  0.00  C4
ATOM     48  C5  XYL      3      4.854 -15.404  -5.899  0.00  0.00  C5
.....
ATOM     74  H1  GLC      4      6.312 -20.297  -3.897  0.00  0.00  H1
ATOM     75  H2  GLC      4      4.081 -21.295  -1.991  0.00  0.00  H2
ATOM     76  H3  GLC      4      4.709 -22.175  -4.882  0.00  0.00  H3
ATOM     77  H4  GLC      4      2.240 -20.654  -3.830  0.00  0.00  H4
ATOM     78  H5  GLC      4      4.504 -19.584  -5.636  0.00  0.00  H5
ATOM     79  H6  GLC      4      1.833 -18.619  -5.118  0.00  0.00  H6
ATOM     80  H24 GLC      4      2.931 -17.520  -4.272  0.00  0.00  H24
.....
ATOM     81  C1  GLC      5      1.057 -21.319  -5.860  0.00  0.00  C1
ATOM     82  C2  GLC      5      0.609 -21.564  -7.318  0.00  0.00  C2
.....
ATOM    199  C1  GAL     11      5.825 -12.737   3.053  0.00  0.00  C1
ATOM    200  C2  GAL     11      4.966 -12.753   1.793  0.00  0.00  C2
ATOM    201  C3  GAL     11      3.486 -12.876   2.114  0.00  0.00  C3
ATOM    202  C4  GAL     11      3.243 -14.031   3.058  0.00  0.00  C4
ATOM    203  C5  GAL     11      4.137 -13.890   4.291  0.00  0.00  C5
ATOM    204  C6  GAL     11      4.003 -15.003   5.307  0.00  0.00  C6
.....
ATOM    255  H5  XYL     13     19.537 -11.651  -0.139  0.00  0.00  H5
ATOM    256  H1  XYL     13     16.585 -11.560  -2.329  0.00  0.00  H1
ATOM    257  H2  XYL     13     18.338 -12.051  -3.974  0.00  0.00  H2
ATOM    258  H3  XYL     13     19.481 -13.841  -1.798  0.00  0.00  H3
ATOM    259  H4  XYL     13     20.596 -11.287  -3.032  0.00  0.00  H4
TER


```

## Color code for monosaccharides

Monosaccharide	Code	Code Couleur
	A2G	
	AAL	
Abequose	ABE	
Aceose	ACE	
Aldotetrose	ALT	
Apicose	API	
Arabinose	ARA	
	BGC	
	BMA	
Acide Dehydroascorbique	DHA	
Fructose	FRU	
Fucose	FUC	
	FUL	
Galactose	GAL	
	GLB	
Glucose	GLC	
Gulose	GUL	
Idose	IDO	
Acide Ketodeoxynonulosonique	KDN	
Acide Ketodeoxyoctulosonique	KDO	
Manose	MAN	
	MMA	
	NAG	
Acide neuraminique N-glycolyl	NEG	
	NDG	
	NGA	
Rhamnose	RHA	
Ribose	RIB	
Acide Sialique	SIA	
Tagatose	TAG	
Talose	TAL	
Xylose	XYL	


## Examples

Particles mode.




Smooth HyperBalls mode


Default Smooth HyperBalls mode




CPK



Wan Der Waals




Licorice




## Texture

CPK




Acid-Basic






**Other examples  
(Renderer Menu)**




Coloration Mode

**GoodSell**







**WaterColor**




**Basic**




**IUPAC**




Chain



Sugar




DNA/RNA




## Secondary Structure Menu


### Secondary Structure




### Ca-Trace




B-factor




Surface Mode




### Surface with static/mobile Cut



### Ring Blending (with sugar coloration)




### SugarRibbons






## Special Effect


**BLUR**




**NOISE**




**BLUR2**




**EDGE**




**VORTEX**



**GRAY**



**GLOW**





## FAQ *(mainly for sugar representation)*

### ***I don't have the surface representation for sugar!***

This is because by default, the surface representation for sugar molecule is disabled. Like that, you can easily represent a glycoprotein in surface mode and locate sugars molecule.

To be able to see the sugar surface, click on **Sugar** toggle button present on the **Surface Menu** (last element of this menu).

### ***I have a sugar molecule, but the SugarRibbons doesn't work!***

It might be several causes for this issue:

1. The 3 letters code for your sugar molecules in your PDB file is not recognize (check *Color code for monosaccharides* page 13).
  - a. Change the 3 letters sugar name in your PDB by a recognized code.
  - b. In **Tune Menu**, uncheck **Sugar Only?** And click on **Enable SugarRibbons** again (maybe twice to uncheck, then check again).
2. You have a connectivity issues, ring of 4 or 7 or more, or 2 ring for 1 sugar.
  - a. No solution here. It's still a beta version and we have to fix that 😊

### **The RingBlending or SugarRibbons coloration is gray!**

The 3 letters code for your sugar molecules in your PDB file is not recognize (check *Color code for monosaccharides* page 13).

### **How to make a screenshot without the UserInterface?**

Press **backspace** touch  (above Return touch).

Then press the **P** touch.

The screenshot is on the same folder of the PDB with the name **testscreen-x.png** with **x** the number of the screenshot. Be careful, when you start a new session of SweetUnityMol, the number **x** will start from 0 again. So your previous screenshot will be replaced (a new version will fix that).

### **I fetch a PDB and when I use the ScreenShot function (touch P), I can't find my screenshot!**

When you fetch a PDB the Screenshot function doesn't work (the screenshot is saved in a folder who doesn't exist). Will be fixed on a next version.

### **I want to have a white background!**

Display ⇒ White

### **I want to use the SugarColorationMode!**

Atom Appareance ⇒ Panels ⇒ Sugar

### **I want to use a Texture on Hyperballs!**

Atom Appareance ⇒ Renderer ⇒ click on a sphere image (you can apply your texture only on a selection with *Residues, Chain, Extended atoms selection*).

**I want to use a Texture on the Surface representation!**

**Surface** ⇒ Generate (if you didn't do it again)

Surface ⇒ Texture ⇒ Click on a sphere image.

**I want to represent a glycoprotein in SurfaceMode and the Sugar molecule in SugarRibbons!**

Surface menu ⇒ Generate

Sugar Menu ⇒ Enable SugarRibbons.

**I want to represent a glycoprotein in Hyperballs for amino acids and SugarRibbons for sugar !**

You have to play with the **Sugar/Non Sugar** checkbox in the **Sugar Menu**. You may have to click twice on these checkboxes (issues with the UserInterface. Will be fixed on a next version).

For example with 3FUS from a new session of SweetUnityMol:

In the input PDB ID : 3FUS ⇒ Fetch PDB

Sugar ⇒ Enable SugarRibbons

Sugar ⇒ Non Sugar (twice)

**Some bonds in my molecule is missing!**

Our shader need OpenGL. Use the file **Run\_SweetUmol\_32.bat** or **Run\_SweetUmol\_32.bat** depending your version (the line has to be something like "umol-winXX-20141002.exe -force-opengl").

**I have a black screen, or a part of the screen stay black**

Do not use the full screen mode (check "**Windowed**" in the "UnityMol Configuration" menu (this menu is the menu where you choose your resolution and click on "Play!").

We discover this issues with the news version of Unity3D. We didn't locate where this problem is from...